


# Tired of dealing with Endotoxin?

## Treat your cells to LoTox™ - Purified Allergens with Low Endotoxin Content!

- Designed for:
  - Human or murine cellular studies
  - T cells, dendritic cells, APC's
  - Histamine Release
  - Animal models of IgE responses
  - Measuring allergen specific IgG and IgE antibodies
- Defined cysteine protease activity (LoTox™ nDer p 1)
- Successfully used in immunotherapy trials (LoTox™ rPhl p 1 and rPhl p 5)

### New LoTox™ Fel d 1 products



- **New LoTox™ products: nDer p 2, nFel d 1, rFel d 1-deglycosylated**

Product Code	Description	Expression System	Unit	Endotoxin*
LTN-DP1-1	LoTox™ nDer p 1	Natural	250µg	<0.03 EU/µg
LTN-DP2-1	LoTox™ nDer p 2	Natural	250µg	<0.03 EU/µg
LTN-FD1-1	LoTox™ nFel d 1	Natural	250µg	<0.03 EU/µg
LTR-FD1-1	LoTox™ rFel d 1	Pichia pastoris	250µg	<0.03 EU/µg
LTR-FD1D-1	LoTox™ rFel d 1-DG (deglycosylated)	Pichia pastoris	250µg	<0.03 EU/µg
LTR-BV1A-1	LoTox™ rBet v 1a	E.coli	250µg	<0.03 EU/µg
LTR-PP1-1	LoTox™ rPh1 p 1	E.coli	250µg	<0.03 EU/µg
LTR-PP5A-1	LoTox™ rPh1 p 5a	E.coli	250µg	<0.03 EU/µg

\* Based on LAL assay

In the U.S.A.: Indoor Biotechnologies Inc. 1216 Harris Street, Charlottesville, VA 22903. Tel: (434) 984-2304; Fax: (434) 984-2709; E-mail: [mail@inbio.com](mailto:mail@inbio.com)

In Europe: Indoor Biotechnologies Ltd. The Old Brewery, 38 High Street, Warminster, Wiltshire, BA12 9AF, UNITED KINGDOM

Tel: 44 (0)5601 153 291 Tel:/FAX: 44 (0)1985 218300; E-mail: [info@indoorbiotech.co.uk](mailto:info@indoorbiotech.co.uk)