

## Natural Mus m 1

**Product Code: NA-MM1-2**

Allergen: Natural Mus m 1

Lot No: xxxxx

Source: Mouse Urine

Mol. Wt: 18-20 kD

Purification: Purified from concentrated, male, mouse urine by gel filtration and ion exchange HPLC. Purity > 95%.

Concentration: See product insert

Formulation: Preservative-free and carrier-free in 20 mM Tris-HCl, pH 8.1. Filtered through 0.22µ filter.

Storage: Store at -20°C

Notes: Not tested for endotoxin



nMus m 1

**For Research Use Only: Not for Diagnostic or Therapeutic Use**

### REFERENCES:

1. Ferrari E, Lodi T, Sorbi RT, Tirindelli R, Cavaggioni A, Spisni A. Expression of a lipocalin in *Pichia pastoris*: secretion, purification and binding activity of a recombinant mouse major urinary protein. *FEBS Lett.* 1997 Jan 13;401(1):73-7.
2. Lorusso, JR, Moffat S, Ohman JL Jr. Immunological and biochemical properties of the major mouse urinary allergen (Mus m 1). *J Allergy Clin Immunol* 1986 Nov;78(5 Pt 1):928-37.