

Natural Rat n 1

Product Code: NA-RN1-2

Allergen:	Natural Rat n 1
Source:	Rat Urine
Mol. Wt:	18-20 kD
Purification:	Purified from concentrated, male, rat urine by gel filtration and ion exchange HPLC. Purity > 95%.
Concentration:	See package insert. By Advanced Protein Assay.
Formulation:	Preservative-free and carrier-free in 20 mM Tris-HCl, pH 8.1. Filtered through 0.22µ filter.
Storage:	Store at -20°C
Notes:	Not tested for endotoxin



For Research Use Only: Not for Diagnostic or Therapeutic Use

REFERENCES:

1. Jeal H, Draper A, Harris J, Taylor AN, Cullinan P, Jones M. Determination of the T cell epitopes of the lipocalin allergen, Rat n 1. Clin Exp Allergy 2004; 34:1919-1925.
2. Bocskei Z, Groom CR, Flower DR, Wright CE, Phillips SE, Cavaggioni A, Findlay JB, North AC. Pheromone binding to two rodent urinary proteins revealed by X-ray crystallography. Nature 1992; 360:186-188.