

## Recombinant Der f 2

**Product Code: RP-DF2-1**

Allergen: Recombinant Der f 2 (*Dermatophagoides farinae* allergen 2)

Lot No: xxxx

Source: *Pichia pastoris*

Mol. Wt: 14kD

Purification: Purified from *Pichia pastoris* culture supernate by mAb affinity chromatography.

Concentration: See product insert

Formulation: Preservative-free and carrier-free in phosphate buffered saline, pH 7.4, sterile filtered.

Storage: Store at -20°C

Notes: (1) The preparation shows a major Der f 2 band at 14kD on silver stained SDS-PAGE having purity > 90%. (2) Not tested for endotoxin.



rDer f 2

**For Research Use Only: Not for Diagnostic or Therapeutic Use**

**REFERENCES:**

1. Heyman PW, Chapman MD, Fox J, Aalberse, RC, PlattsMills, TAE. Antigenic and structural analysis of the Group II allergens (Derp II and Der f II) from house dust mites (*Dermatophagoides* sp). *J Allergy Clin Immunology* 1989;83:1055-67.
2. Mueller GA, Smith AM, Williams DC Jr., Hakkaart GAJ, Aalberse RC, Chapman MD, Rule, GS, Benjamin DC. Expression and secondary structure determination by NMR methods of the major house dust mite allergen Der p 2. *J Biol Chem* 1997;272:26893-8.
3. Smith AM, Benjamin DC, Hozic N, Derewenda U, Smith WA, Thomas WR, Gafvelin G, van Hage-Hamsten M, Chapman MD. The molecular basis of antigenic cross-reactivity between the group 2 mite allergens. *J Allergy Clin Immunol* 2001;107:977-84.